




Product Datasheet

GTF2I antibody (orb37140)

Description	Rabbit polyclonal antibody to GTF2I is a multifunctional phosphoprotein with roles in transcription ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	GTF2I
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	112416
Uniprot ID	P78347
NCBI	NP_001267729.1 , NP_127492.1 , NP_127493.1 , NP_001509.3 , NP_001157108.1 , NP_127494.1
Dilution Range	IF/ICC: 1:10-50, WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of K562 cells (...)



Western blot analysis of Hela cell line ...