

Product Datasheet

GALE antibody (orb37163)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, IHC-P, WB

Immunogen Synthetic Peptide

Target GALE

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 38282

Uniprot ID [Q14376](#)


NCBI [NP_001008217.1](#), [NP_001121093.1](#), [NP_000394.2](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50

Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of WiDr cell...



Immunohistochemical staining of paraffin...