


Product Datasheet

TCP1 theta/CCT8 Antibody (orb371727)


Description	TCP1 theta/CCT8 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IHC-Fr, WB
Immunogen	E.coli-derived human CCT8 recombinant protein (Position: K421-D548). Human CCT8 shares 98.4% amino acid (aa) sequence identity with mouse CCT8.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59621 MW
Uniprot ID	P50990
Expiration Date	12 months from date of receipt.



WB analysis of CCT8 using anti-CCT8 anti...



IHC analysis of CCT8 using anti-CCT8 ant...



IHC analysis of CCT8 using anti-CCT8 ant...