

Product Datasheet

BCA1/CXCL13 Antibody (orb371733)

Description	BCA1/CXCL13 Antibody	
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	Western blot analysis of recombinant mouse...
Immunogen	E. coli-derived mouse BCA1 recombinant protein (Position: I22-A109). Mouse BCA1 shares 45.3% amino acid (aa) sequence identity with human BCA1.	
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	Immunohistochemical staining of mouse ly...
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	11927 MW	
Uniprot ID	O55038	
Expiration Date	12 months from date of receipt.	