



## Product Datasheet

### DDAH1 Antibody (orb371734)


|                            |   |
|----------------------------|---|
| <b>Description</b>         | DDAH1 Antibody  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | FC, ICC, IF, IHC, WB  |
| <b>Immunogen</b>           | A synthetic peptide corresponding to a sequence at the C-terminus of human DDAH1 (195-226aa QKALKIMQMSDHRVDKLTVPDDIAANCIYLN), different from the related mouse and rat sequences by one amino acid.   |
| <b>Form/Appearance</b>     | Lyophilized   |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.   |
| <b>Storage</b>             | Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Add 0.2ml of distilled water will yield a concentration of 500ug/ml.  |
| <b>Isotype</b>             | Rabbit IgG  |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 31122 MW  |
| <b>Uniprot ID</b>          | <a href="#">O94760</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Flow Cytometry analysis of A431 cells us...



WB analysis of DDAH1 using anti-DDAH1 an...



IF analysis of DDAH1 using anti-DDAH1 an...