

Product Datasheet

Mouse IL1 Receptor protein (orb371811)

Description	Recombinant Mouse Interleukin-1 receptor type 1 is produced by our Mammalian expression system and t...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Application notes	(C-Fc)
Uniprot ID	P13504
Expiration Date	6 months from date of receipt.