

Product Datasheet

Human Cytotoxic T-lymphocyte 4 protein (orb371903)

Description	Recombinant Human Cytotoxic T-lymphocyte Protein 4 is produced by our Mammalian expression system an...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Flag
Note	For research use only
Application notes	(C-Flag)
Uniprot ID	P16410
Expiration Date	6 months from date of receipt.