

Product Datasheet

Human S100A13 protein (orb371909)

Description	Recombinant Human Protein S100-A13 is produced by our E.coli expression system and the target gene e...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mMTris HCl,1mM CaCl_2 ,0.1%Tween-20,pH8.0.
Storage	Lyophilized protein should be stored at -20°C , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Uniprot ID	Q99584
Expiration Date	6 months from date of receipt.