


Product Datasheet

NLRP1 antibody (orb37218)

Description	NLRP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 1 protein (1-146AA)
Target	NLRP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9C000
Dilution Range	IF:1:100-200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of SH-SY5Y C...