

Product Datasheet

MMP-3 Inhibitor I (orb372565)

Description	Peptide to MMP-3 Inhibitor I
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Arg-Cys-Gly-Val-Pro-Asp-NH ₂
Formula	C ₂₇ H ₄₆ N ₁₀ O ₉ S
Purity	> 95%
Source	Synthetic
MW	686.80
CAS Number	158841-76-0
Expiration Date	12 months from date of receipt.

