

Product Datasheet

Papain Inhibitor (orb372660)

Description	Peptide to Papain Inhibitor
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Gly-Tyr-Arg
Formula	C ₁₉ H ₂₉ N ₇ O ₆
Purity	> 95%
Source	Synthetic
MW	451.49
CAS Number	70195-20-9
Expiration Date	12 months from date of receipt.

