

Product Datasheet

TRAP-6 (2-6) (orb372881)

| | |
|-------------------------|--|
| Description | Peptide to TRAP-6 (2-6) |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Protein Sequence | Phe-Leu-Leu-Arg-Asn |
| Formula | $C_{31}H_{51}N_9O_7$ |
| Purity | > 95% |
| Source | Synthetic |
| MW | 661.81 |
| CAS Number | 141136-84-7 |
| Expiration Date | 12 months from date of receipt. |

