

Product Datasheet

Triptorelin (orb372886)

| | |
|-------------------------|---|
| Description | Peptide to Triptorelin |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Protein Sequence | Glp-His-Trp-Ser-Tyr-D-Trp-Leu-Arg-Pro-Gly-NH ₂ |
| Formula | C ₆₄ H ₈₃ N ₁₈ O ₁₃ |
| Purity | > 95% |
| Source | Synthetic |
| MW | 1311.50 |
| CAS Number | 57773-63-4 |
| Expiration Date | 12 months from date of receipt. |

