


Product Datasheet

ERVFRD-1 antibody (orb373072)

Description	Rabbit polyclonal antibody to ERVFRD-1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ERVFRD-1(NP_997465.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60KDa
Uniprot ID	P60508
Dilution Range	WB: 1:1000-4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...