



Product Datasheet

IL20 antibody (orb373436)


Description	Rabbit polyclonal antibody to IL20
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human IL20(NP_061194.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20kDa
Uniprot ID	Q9NYY1
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human li...



Immunohistochemical staining of mouse lu...



Immunohistochemical staining of mouse he...