



Product Datasheet

AGO2 antibody (orb373479)


Description	Rabbit polyclonal antibody to AGO2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Argonaute-2(NP_036286.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110kDa
Uniprot ID	Q9UKV8
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-100, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Immunohistochemical staining of mouse li...



Western blot analysis of extracts of var...