

Product Datasheet

GLUD1 antibody (orb373657)

Description Rabbit polyclonal antibody to GLUD1

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, IP, WB

Immunogen Recombinant fusion protein of human GLUD1(NP_005262.1).

Preservatives PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 52kDa


Uniprot ID **P00367**

Dilution Range WB: 1:500-2000, IHC-P: 1:100-200


Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of human li...



Immunohistochemical staining of human co...