



Product Datasheet

NOTCH1 antibody (orb373660)

Description	Rabbit polyclonal antibody to NOTCH1, which is also known as AOVD. NOTCH1 is frequently found in the...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Notch1(NP_060087.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	300KDa
Uniprot ID	P46531
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human esophagus tissue.



Western blot analysis of extracts of HT-29 cells.