

Product Datasheet

TRPC3 antibody (orb373732)

Description Rabbit polyclonal antibody to TRPC3

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen A synthetic peptide of human TRPC3(NP_001124170.1).

Preservatives PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 96kDa


Uniprot ID [Q13507](#)

Dilution Range WB: 1:500-2000, IHC-P: 1:50-200


Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of mouse br...



Immunohistochemical staining of rat brai...