



Product Datasheet

ARHGEF11 antibody (orb373879)


Description	Rabbit polyclonal antibody to ARHGEF11
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human ARHGEF11
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	168kDa
Uniprot ID	O15085
Entrez	9826
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...



Immunohistochemical staining of human pr...



Immunohistochemical staining of human lu...