



Product Datasheet

MYLK antibody (orb373925)


Description	Rabbit polyclonal antibody to MYLK
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human MYLK(NP_444253.3).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16kDa/80kDa/110kDa/197-210kDa
Uniprot ID	Q15746
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Mou...



Immunohistochemical staining of human ga...



Immunohistochemical staining of rat kidn...