



Product Datasheet

CETN3 antibody (orb373959)


| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to CETN3 |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, IHC, WB |
| Immunogen | Recombinant fusion protein of human CETN3(NP_004356.2). |
| Preservatives | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 20kDa |
| Uniprot ID | O15182 |
| Dilution Range | WB: 1:500-2000, IHC-P: 1:50-100 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extract of vari...



Immunohistochemical staining of human li...



Immunohistochemical staining of human lu...