

Product Datasheet

PTTG1 antibody (orb374037)

Description Rabbit polyclonal antibody to PTTG1

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Recombinant fusion protein of human PTTG1(NP_001269311.1).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 24kDa


Uniprot ID [O95997](#)

Dilution Range WB: 1:500-2000, IHC-P: 1:100-200


Expiration Date 12 months from date of receipt.



Western blot analysis of extract of vari...



Immunohistochemical staining of human lu...



Immunohistochemical staining of human li...