



Product Datasheet

WNT5B antibody (orb374043)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to WNT5B |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | Recombinant fusion protein of human WNT5B(NP_116031.1). |
| Preservatives | PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 50kDa |
| Uniprot ID | Q9H1J7 |
| Dilution Range | WB: 1:500-2000, IHC-P: 1:50-100 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extract of vari...



Immunohistochemical staining of mouse lu...