

Product Datasheet

GNE antibody (orb374202)

Description Rabbit polyclonal antibody to GNE

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Recombinant fusion protein of human GNE(NP_001121699.1).

Preservatives PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 79kDa


Uniprot ID [Q9Y223](#)

Dilution Range WB: 1:500-2000, IHC-P: 1:50-100


Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of mouse br...



Immunohistochemical staining of rat brai...