


Product Datasheet


ATP1B2 antibody (orb374297)

Description	Rabbit polyclonal antibody to ATP1B2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 68-290 of human ATP1B2 (NP_001669.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42kDa
Uniprot ID	P14415
Dilution Range	WB: 1:200-500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells

...



Western blot analysis of extracts of var...