


Product Datasheet

PRKACA (Phospho-T197) antibody (orb374347)

Description	Rabbit polyclonal antibody to PRKACA (Phospho-T197)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T197 of human PKA C-alpha (PRKACA)(NP_002721.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa
Uniprot ID	P17612
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of NIH...