


Product Datasheet

PDGFRA (Phospho-Y849) antibody (orb374348)

Description	Rabbit polyclonal antibody to PDGFRA (Phospho-Y849)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y849 of human PDGFRA(NP_006197.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	190kDa
Uniprot ID	P16234
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Hep...