

Product Datasheet

PKC (pan) (gamma) (Phospho-T514) antibody (orb374354)

Description	Rabbit polyclonal antibody to PKC (pan) (gamma) (Phospho-T514)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T514 of human PKC (pan) (gamma)PRKCA(NP_002728.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80kDa
Uniprot ID	P05129 , P17252 , Q02156 , Q05655 , P24723 , P05771 , Q04759
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.

