



Product Datasheet

PCSK9 antibody (orb374429)


Description	Rabbit polyclonal antibody to PCSK9
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human Proprotein Convertase 9(PCSK9)(NP_777596.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80KDa
Uniprot ID	Q8NBP7
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse li...



Immunohistochemical staining of human st...



Immunohistochemical staining of rat kidn...