



Product Datasheet

CACYBP antibody (orb374468)


Description	Rabbit polyclonal antibody to CACYBP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CACYBP(NP_055227.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	Q9HB71
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of NIH-3T3 C...



Western blot analysis of extracts of var...