



## Product Datasheet

### PMAIP1 antibody (orb374505)


<b>Description</b>	Rabbit polyclonal antibody to PMAIP1
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IF, WB
<b>Immunogen</b>	A synthetic peptide of human PMAIP1(NP_066950.1).
<b>Preservatives</b>	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	14kDa
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/Q13794">Q13794</a>
<b>Dilution Range</b>	WB: 1:500-2000
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of mouse bre...



Immunofluorescence analysis of mouse bre...



Immunofluorescence analysis of rat breas...