

Product Datasheet

MTHFD2 antibody (orb374613)

Description Rabbit polyclonal antibody to MTHFD2

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen Recombinant fusion protein of human MTHFD2(NP_006627.2).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 38KDa


Uniprot ID **P13995**

Dilution Range WB: 1:500-2000


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of mouse kidney



Immunohistochemical staining of rat kidney



Immunohistochemical staining of human liver