

## Product Datasheet

### CYP11A1 peptide (orb374896)

<b>Description</b>	CYP11A1 peptide
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HPLC
<b>Target</b>	CYP11A1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb313866 (CYP11A1 antibody)
<b>Source</b>	Synthesis
<b>Expiration Date</b>	6 months from date of receipt.