
Product Datasheet

**Adipophilin Antibody / Perilipin 2 / ADRP / PLIN2
(orb376629)**

Description

Adipophilin, also called Perilipin 2, ADRP (Adipose

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FACS, WB

Immunogen

A human recombinant partial protein was used as the immunogen for this Adipophilin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

The Adipophilin antibody (with azide) can be stored at 2-8°C. The azide-free format should be aliquoted and stored at -20°C or colder.

Note

For research use only

Application notes

Titering of the Adipophilin antibody may be required for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID
[Q99541](#)
Hazard Information


This Adipophilin antibody is available for research use only.

Entrez
[123](#)
Dilution Range


Flow cytometry: 0.5-1ug/million cells,Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.



Western blot testing of 1) partial recom...



FACS testing of human PBMC with Adipophi...