
Product Datasheet

TNFSF15 Antibody / VEGI / TL1A (orb376640)

Description

Vascular endothelial growth inhibitor (VEGI), also


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

IHC-P



IHC testing of FFPE human parathyroid tu...

Immunogen

A full length human recombinant protein was used as the immunogen for this TL1A antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage


Store the TL1A antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the TL1A antibody to be titrated up or down for optimal performance.



IHC testing of FFPE human colon carcinom...

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[O95150](#)

Hazard Information

This TL1A antibody is available for research use only.

Entrez

[9966](#)

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.