

Product Datasheet

NIPSNAP3B antibody (FITC) (orb37730)

Description	NIPSNAP3B antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Protein NipSnap homolog 3B protein (1-247AA)
Target	NIPSNAP3B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BS92
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.