



Product Datasheet

AKAP13 antibody (orb377820)

Description	Mouse monoclonal antibody to AKAP13.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	AKAP13
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	307550
Uniprot ID	Q12802
Dilution Range	WB: 1:8000
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...