

Product Datasheet

LIAS antibody (orb377861)

Description	Rabbit polyclonal antibody to LIAS.
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	This LIAS antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 298-330 amino acids from human LIAS.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	O43766 , Q99M04 , Q5XIH4
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

