


Product Datasheet

ACKR1 antibody (Biotin) (orb377882)

Description	Goat polyclonal antibody against ACKR1 conjugated to Biotin.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Target	Duffy / FY / DARC, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 37kDa band observed in Human Liver lysates (calculated MW of 35.7kDa according to NP_001116423.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	HRAELSPSTENSSQLDFEDC.
Clonality	Polyclonal
MW	35.7; 35.6
Entrez	2532
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...