




Product Datasheet

PCSK9 antibody (orb37796)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PCSK9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74286
Uniprot ID	Q8NBP7
NCBI	NP_777596.2
Dilution Range	FACS: 1:10-50, IHC-P: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2, HeLa, SW...



Flow cytometric analysis of HeLa cells u...