

Product Datasheet

CD20 antibody (orb378213)

Description	Recombinant monoclonal antibody to CD20
Species/Host	Human
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human lymphoblastoid cell line SB.
Target	CD20
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P11836
Expiration Date	12 months from date of receipt.