

Product Datasheet

OX40L antibody (orb378214)

Description	Recombinant monoclonal antibody to OX40L
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	Oxelumab binds to human OX40L, the ligand of OX40, stimulatory receptor that is expressed on activated T cells, natural killer (NK) cells and natural killer T (NKT) cells.
Target	OX40L
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to the co-stimulatory OX40L, inhibiting its interaction with OX40.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P23510
Expiration Date	12 months from date of receipt.