

Product Datasheet

CD52 antibody (orb378215)

Description	Recombinant monoclonal antibody to CD52
Species/Host	Rat
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, WB
Immunogen	Human peripheral mononuclear cells enriched for T-cells.
Target	CD52
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P31358
Expiration Date	12 months from date of receipt.