


## Product Datasheet

### Sp3/4 antibody (orb378416)

<b>Description</b>	Rabbit polyclonal antibody to Sp3/4
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB: 1:1000 IHC-P: 1:200
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...