




Product Datasheet

GIT1-Y554 antibody (orb37916)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	GIT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	84341
Uniprot ID	Q9Y2X7
NCBI	NP_001078923.1 , NP_054749.2
Dilution Range	WB 1:50-1:100
Expiration Date	12 months from date of receipt.



WB analysis of mouse NIH-3T3 cell line l...



WB analysis of U-2OS whole cell lysate (...)