



Product Datasheet

BNIP3 Antibody (orb379656)


| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to BNIP3 |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC, ICC, IHC |
| Immunogen | Purified recombinant fragment of human BNIP3 (AA: 50-155) expressed in E. Coli. |
| Preservatives | Purified antibody in PBS with 0.05% sodium azide |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | Mouse IgG2a |
| Clonality | Monoclonal |
| MW | 21.5kDa |
| Entrez | 664 |
| Dilution Range | IHC: 1/200 - 1/1000; ICC: 1/50 - 1/250; RCM: 1/200 - 1/400; ELISA: 1/10000 |
| Expiration Date | 12 months from date of receipt. |



Line graph illustrates about the Ag-Ab r...



Flow cytometric analysis of HeLa cells u...



Immunofluorescence analysis of HeLa cell...