



Product Datasheet

UCP3 Antibody (orb379666)


Description	Mouse monoclonal antibody to UCP3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Purified recombinant fragment of human UCP3 (AA: 1-113 and 217-312) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	34.2kDa
Entrez	7352
Dilution Range	ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Immunofluorescence analysis of HL-7702 C...



Western blot analysis of HEK293 cell lys...