

## Product Datasheet

### GFP-Tag Monoclonal Antibody(1G6) (orb380378)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | GFP-Tag Monoclonal Antibody(1G6)  |
| <b>Species/Host</b>        | Mouse   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Monoclonal  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |