

Product Datasheet

HSPB8/HSP22 Monoclonal Antibody(2C3) (orb380460)

Description	HSPB8/HSP22 Monoclonal Antibody(2C3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q9UJY1
Expiration Date	12 months from date of receipt.